

PRIDE
MOBILITY®

TRAVEL POWER CHAIR

EXPLORE THE WORLD

Whether you are riding the rails or cruising the seas, you need a travel power wheelchair that is lightweight and portable. Consider a Pride travel power wheelchair for your next adventure.



EASY DISASSEMBLY
INTO 5 LIGHTWEIGHT PIECES



EASY FOLD DESIGN
FOR PORTABILITY AND STORAGE



SCAN THE QR CODE
TO LEARN MORE ABOUT
TRAVEL MOBILITY.

   
pridemobility.com

PRIDE
MOBILITY®

GO CHAIR®



FDA Class II
Medical Device¹

FEATURES

- Easily disassembles into 5 lightweight pieces
- Dual swing-away, under seat storage bins
- Number one for portability
- Heaviest piece is only 36 lbs.
- Available in 6 colors

ACCESSORIES

- Armrest mount cane clip
- Cane/crutch holder (3" single, 3" double)
- Cup holder
- Lap belt (50", 70")
- Oxygen tank holder
- Rear mount bracket
- Saddlebag
- Safety flag
- Weather cover
- XLR USB charger
- Additional battery pack



COLORS

Ruby Red	Pearl White
Amber Orange	Citrine Yellow
Rose Quartz	Sapphire Blue

MILITARY PATCHES



300lbs.
Weight Capacity

Up To
3.8mph
Maximum Speed

25.5"
Turning Radius

Up To
13.2
Miles Per Charge²

JAZZY® ULTRA LIGHT



FDA Class II
Medical Device¹

FEATURES

- Easy fold design that locks into a folded position
- Weighs just 33 pounds without the battery
- Mesh back and seat with standard cushion for extra support
- Programmable joystick controller can attach to either armrest
- Powered by an airline compliant³ lithium battery (10Ah Only)

ACCESSORIES

- Jazzy® Travel Bag
- Additional battery
- Headrest & attendant control bracket bundle
- Padded armrest covers



COLORS

Black

FOLDED SIZE
24.5" x 31.5" x 14"



300lbs.
Weight Capacity

Up To
3.7mph
Maximum Speed

33.19"
Turning Radius

Up To
16.4
Miles Per Charge²

WARRANTY

FRAME: 5-year limited
ELECTRONICS: 13-month limited
DRIVETRAIN: 13-month limited
BATTERY: 13-month limited

HCPCS Code: K0899

1. FDA Class II Medical Devices designed to aid individuals with mobility impairments.
2. Range varies based on user weight, terrain and battery type.
3. Lithium-ion battery is subject to federal regulations and airline policy. Please consult your airline prior to flight.



pridemobility.com

**PRIDE
MOBILITY®**

AUTHORIZED PROVIDER